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VR	Name	Beschreibung
IS	Integer String	A string of characters representing an Integer in base-10 (decimal), shall contain only the characters 0 - 9, with an optional leading "+" or "-". It may be padded with leading and/or trailing spaces. Embedded spaces are not allowed. The integer, n, represented shall be in the range: $-2^{31} \leq n \leq (2^{31} - 1)$. 12 Bytes maximum.
LO	Long String	A character string that may be padded with leading and/or trailing spaces. The character code 5CH (the BACKSLASH) shall not be present, as it is used as the delimiter between values in multiple valued Default Character Repertoire and/or as defined by (0008,0005). 64 chars maximum.
PN	Person Name	A character string encoded using a 5 component convention. The character code 5CH (the BACKSLASH) shall not be present, as it is used as the delimiter between values in multiple valued data elements. The string may be padded with trailing spaces. The five components in their order of occurrence are: family name complex, given name complex, middle name, name prefix, name suffix. Any of the five components may be an empty string. The component delimiter shall be the caret '^' character (5EH). Delimiters are required for interior null components. Trailing null components and their delimiters may be omitted. Multiple entries are permitted in each component and are encoded as natural text strings, in the format preferred by the named person. This conforms to the ANSI HISPP MSDS Person Name common data type. 64 character maximum per component group.
SH	Short String	A character string that may be padded with leading and/or trailing spaces. The character code 05CH (the BACKSLASH) shall not be present, as it is used as the delimiter between values for multiple data elements. The string

		string shall not have Control Characters except ESC. Default Character Repertoire and/or as defined by (0008,0005). 16 chars maximum.
UI	Unique Identifier (UID)	A character string containing a UID that is used to uniquely identify a wide variety of items. The UID is a series of numeric components separated by the period "." character. If a Value Field containing one or more UIDs is an odd number of bytes in length, the Value Field shall be padded with a single trailing NULL (00H) character to ensure that the Value Field is an even number of bytes in length. 64 chars maximum.

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Usage:

dcm2pgm [commands] <source>

Converts a DICOM file to a PGM or PPM file, depending on the contents of the source. By default only default metadata required for a Secondary Capture image are written besides the image data.

<source> the file or directory containing the DICOM file(s).

Commands:

-a <study_id> Set the StudyID tag (00200010) to <study_id>.
-A <study_uid> Set the StudyInstanceUID tag (0020000d) to <study_uid>.
-b <series_num> Set the SeriesNumber tag (00200011) to <series_num>.
-B <series_uid> Set the SeriesInstanceUID tag (0020000e) to <series_uid>.
-c <image_num> Set the InstanceNumber tag (00200013) to <image_num>.
-C <sop_uid> Set the SOPInstanceUID tag (00080018) to <sop_uid>.
-d <file-uri> The name of the destination file or a directory.
-F Include all DICOM tags as image comments.
-h Print this help message.
-n <name> Set the PatientName tag (00100010) to <name>.
-N <pat_id> Set the PatientID tag (00100020) to <pat_id>.
-x Include only DICOM tags required for a Secondary Capture image.
 If not enough information is given in the base image, default
 values are added.
-X Include no DICOM tags in the image.

<file-uri> Describes a file in a operating-system independent way. See the
 API-Doc of the URI class. For Windows-OS the absolute URI
 file:/c:/user/tom/foo.txt
 describes the file
 C:\user\tom\foo.txt
 Relative URI's, e.g. without the 'file:' schema-prefix, are
 relative to the user-directory, given by the system property
 user.dir. For example: If the user.dir is
 C:\user\tom\
 and the relative URI is
 /abc/foo.txt
 the referenced file is
 C:\user\tom\abc\foo.txt
 The abbreviations '.' for the current and '..' for the upper
 directory are valid to form a relative URI.

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This library is free software; you can redistribute it and/or modify it under the terms of the GNU General Public License as published by the Free Software Foundation; either version 2 of the License, or (at your option) any later version.

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pgm2dcm version 1.0

Usage:
pgm2dcm [commands] <source>

Converts a PGM or PPM file to DICOM. By default only default metadata required for a Secondary Capture image are written besides the image data.

<source> the file or directory containing the DICOM file(s).

Commands:

-a <study_id> Set the StudyID tag (00200010) to <study_id>.
-A <study_uid> Set the StudyInstanceUID tag (0020000d) to <study_uid>.
-b <series_num> Set the SeriesNumber tag (00200011) to <series_num>.
-B <series_uid> Set the SeriesInstanceUID tag (0020000e) to <series_uid>.
-c <image_num> Set the InstanceNumber tag (00200013) to <image_num>.
-C <sop_uid> Set the SOPInstanceUID tag (00080018) to <sop_uid>.
-d <file-uri> The name of the destination file or a directory.
-f Write all DICOM tags found in the base image. Add default metadata required for a Secondary Capture image. Preserve the following tags:
 PatientName (00100010)
 PatientID (00100020)
 PatientBirthDate (00100030)
 PatientSex (00100040)
-F Write only DICOM tags found in the base image.
-h Print this help message.
-n <name> Set the PatientName tag (00100010) to <name>.
-N <pat_id> Set the PatientID tag (00100020) to <pat_id>.

<file-uri> Describes a file in a operating-system independend way. See the API-Doc of the URI class. For Windows-OS the absolute URI
 file:/c:/user/tom/foo.txt
 describes the file
 C:\user\tom\foo.txt
Relative URI's, e.g. without the 'file:' schema-prefix, are
relative to the user-directory, given by the system property
user.dir. For example: If the user.dir is
 C:\user\tom\
and the relative URI is
 /abc/foo.txt
the referenced file is
 C:\user\tom\abc\foo.txt
The abbreviations '.' for the current and '..' for the upper
directory are valid to form a relative URI.

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IFTM Institut fuer Telematik in der Medizin GmbH, www.iftm.de